

# Contents

|   |          |
|---|----------|
| <b>OH Parser - Complete Documentation</b>                 | <b>1</b> |
| Table of Contents . . . . .                               | 1        |
| Overview . . . . .  | 2        |
| Key Features . . . . .                                    | 2        |
| Design Principles . . . . .                               | 2        |
| Project Structure . . . . .                               | 2        |
| Installation . . . . .                                    | 2        |
| Quick Start . . . . .                                     | 3        |
| OH Profile JSON Structure . . . . .                       | 3        |
| Top-Level Structure . . . . .                             | 3        |
| EMG Nested Structure (3 Levels) . . . . .                 | 4        |
| API Reference . . . . .                                   | 4        |
| Loading Functions . . . . .                               | 4        |
| Inspection Functions . . . . .                            | 4        |
| Extraction Functions . . . . .                            | 5        |
| Path Resolution Functions . . . . .                       | 6        |
| Filtering . . . . .                                       | 6        |
| Usage Examples . . . . .                                  | 7        |
| Example 1: Weekly EMG Summary (Wide Format) . . . . .     | 7        |
| Example 2: All Session-Level Data (Long Format) . . . . . | 7        |
| Example 3: Filter by Subjects and Date Range . . . . .    | 7        |
| Example 4: Check Data Availability . . . . .              | 8        |
| Example 5: Manual Path Navigation . . . . .               | 8        |
| Path Syntax & Wildcards . . . . .                         | 8        |
| Dot Notation . . . . .                                    | 8        |
| Wildcards in <code>value_paths</code> . . . . .           | 8        |
| Level Names (Implicit Wildcards) . . . . .                | 9        |
| Exclusion Patterns . . . . .                              | 9        |
| CLI Interface . . . . .                                   | 9        |
| EMG Data Reference . . . . .                              | 9        |
| Metric Groups . . . . .                                   | 9        |
| Aggregation Levels . . . . .                              | 10       |
| Pipeline Architecture . . . . .                           | 10       |
| Data Flow . . . . .                                       | 11       |
| Notes . . . . .   | 11       |

## OH Parser - Complete Documentation

Comprehensive guide for extracting data from Occupational Health (OH) profile JSON files into pandas DataFrames.

---

### Table of Contents

1. Overview
2. Project Structure
3. Installation
4. Quick Start
5. OH Profile JSON Structure
6. API Reference
7. Usage Examples

8. Path Syntax & Wildcards
  9. CLI Interface
  10. EMG Data Reference
  11. Pipeline Architecture
- 

## Overview

**OH Parser** is a Python library for extracting data from Occupational Health (OH) profile JSON files into pandas DataFrames for statistical analysis.

## Key Features

- **Load** JSON profiles from a directory (`load_profiles`)
- **Inspect** profile structure with tree visualization (`inspect_profile`)
- **Extract** data using dot-notation paths:
  - `extract()` -> Wide format (one row per subject)
  - `extract_nested()` -> Long format (one row per session/date/side)
  - `extract_flat()` -> Flatten all nested keys
- **Filter** by subjects, date ranges, groups, or data availability (`create_filters`)
- **CLI** for quick exploration without writing code

## Design Principles

1. **Pure function-based** - No classes, just functions and dictionaries
  2. **Dot-notation paths** - Navigate nested JSON with strings like "a.b.c"
  3. **Wildcard support** - Extract all keys with "EMG\_intensity.\*"
  4. **Minimal dependencies** - Only requires pandas
- 

## Project Structure

```
oh_parser/  
|-- __init__.py      # Public API exports  
|-- __main__.py      # CLI entrypoint (python -m oh_parser)  
|-- cli.py           # Command-line interface  
|-- loader.py        # Load OH profile JSON files  
|-- path_resolver.py # Dot-notation path navigation  
|-- filters.py       # Subject/date filtering  
|-- extract.py       # DataFrame extraction functions  
+-- utils.py         # Utility functions
```

---

## Installation

Copy the `oh_parser/` folder to your project, then install dependencies:

```
pip install pandas
```

Or with requirements file:

```
pip install -r requirements.txt
```

**Dependencies:** pandas >= 1.5.0, Python >= 3.9

---

## Quick Start

```
from oh_parser import load_profiles, list_subjects, inspect_profile, extract, extract_nested

# 1. Load all profiles
profiles = load_profiles("/path/to/OH_profiles/")

# 2. List subjects
subjects = list_subjects(profiles) # ['80', '81', '82', ...]

# 3. Inspect structure
inspect_profile(profiles[subjects[0]], max_depth=4)

# 4. Extract specific values (wide format - one row per subject)
df = extract(profiles, paths={
    "p50_L": "sensor_metrics.emg.EMG_weekly_metrics.left.EMG_apdf.active.p50",
    "p50_R": "sensor_metrics.emg.EMG_weekly_metrics.right.EMG_apdf.active.p50",
})

# 5. Extract nested data (long format - one row per session)
df = extract_nested(
    profiles,
    base_path="sensor_metrics.emg",
    level_names=["date", "session", "side"],
    value_paths=["EMG_intensity.*", "EMG_rest_recovery.*"],
    exclude_patterns=["EMG_daily_metrics", "EMG_weekly_metrics"],
)
```

---

## OH Profile JSON Structure

Each subject has one file: {subject\_id}\_OH\_profile.json

### Top-Level Structure

```
{
  "meta_data": { "age": 45, "gender": "M", ... },
  "single_instance_questionnaires": {
    "personal": { ... },
    "biomechanical": { ... },
    "psychosocial": { ... },
    "environmental": { ... }
  },
  "daily_questionnaires": {
    "workload": { "YYYY-MM-DD": { ... } },
    "pain": { "YYYY-MM-DD": { ... } }
  },
  "sensor_metrics": {
    "sensor_timeline": { ... },
    "human_activities": { ... },
    "heart_rate": { ... },
    "posture": { ... },
    "noise": { ... },
    "emg": { ... },
  }
}
```

```

    "wrist_activities": { ... }
  }
}

```

## EMG Nested Structure (3 Levels)

```

sensor_metrics.emg/
|-- {date: DD-MM-YYYY}/          # Level 1: Recording day
|  |-- {session: HH-MM-SS}/      # Level 2: Session start time
|  |  |-- left/                 # Level 3: Side
|  |  |  |-- EMG_session/       # Session metadata
|  |  |  |-- EMG_intensity/     # Intensity metrics
|  |  |  |-- EMG_apdf/         # APDF percentiles
|  |  |  |-- EMG_rest_recovery/ # Rest/recovery metrics
|  |  |  +-- EMG_relative_bins/ # Relative intensity bins
|  |  +-- right/
|  |  +-- ... (same structure)
|  +-- EMG_daily_metrics/       # Aggregated daily
|      |-- left/ { ... }
|      +-- right/ { ... }
+-- EMG_weekly_metrics/        # Aggregated weekly
    |-- left/ { ... }
    +-- right/ { ... }

```

---

## API Reference

### Loading Functions

**load\_profiles(directory, subject\_ids=None, verbose=True) -> Dict[str, dict]** Load all OH profiles from a directory.

| Parameter   | Type        | Description                              |
|-------------|-------------|--|
| directory   | str \  Path | Path to directory containing OH profiles |
| subject_ids | List[str]   | Optional: only load specific subjects    |
| verbose     | bool        | Print loading progress                   |

```

profiles = load_profiles("/path/to/OH_profiles/")
profiles = load_profiles("/path/to/OH_profiles/", subject_ids=["103", "104"])

```

**list\_subjects(profiles) -> List[str]** Get sorted list of subject IDs (sorted numerically).

**get\_profile(profiles, subject\_id) -> dict | None** Get a single profile by subject ID. Returns None if not found.

**load\_profile(filepath) -> dict** Load a single OH profile JSON file.

### Inspection Functions

**inspect\_profile(profile, base\_path="", max\_depth=4, show\_values=False)** Pretty-print the structure of a profile as a tree.

```
inspect_profile(profile)
inspect_profile(profile, base_path="sensor_metrics.emg", max_depth=3)

get_available_paths(profile, base_path="", max_depth=6) -> List[str] Get all dot-notation
paths available in a profile.

paths = get_available_paths(profile)
# ['meta_data.age', 'sensor_metrics.emg.EMG_weekly_metrics.left.EMG_apdf.active.p50', ...]

summarize_profiles(profiles) -> pd.DataFrame Generate summary showing data availability across
subjects.

summary = summarize_profiles(profiles)
# Columns: subject_id, has_meta_data, has_emg, has_EMG_weekly_metrics, ...
```

## Extraction Functions

**extract(profiles, paths, filters=None, include\_subject\_id=True)** -> **pd.DataFrame** Extract specific paths into **wide-format DataFrame** (one row per subject).

| Parameter          | Type            | Description   |
|--------------------|-----------------|---|
| profiles           | Dict[str, dict] | Loaded profiles                                     |
| paths              | Dict[str, str]  | Mapping: column name -> dot-notation path           |
| filters            | dict            | Optional filters from <code>create_filters()</code> |
| include_subject_id | bool            | Include <code>subject_id</code> column              |

```
df = extract(profiles, paths={
    "age": "meta_data.age",
    "emg_p50": "sensor_metrics.emg.EMG_weekly_metrics.left.EMG_apdf.active.p50",
})
# Output: subject_id | age | emg_p50
```

**extract\_nested(profiles, base\_path, level\_names, value\_paths=None, filters=None, exclude\_patterns=None, flatten\_values=True)** -> **pd.DataFrame** Extract nested structures into **long-format DataFrame** (one row per leaf node).

| Parameter        | Type            | Description   |
|------------------|-----------------|---|
| profiles         | Dict[str, dict] | Loaded profiles                                     |
| base_path        | str             | Starting path (e.g., "sensor_metrics.emg")          |
| level_names      | List[str]       | Names for nesting levels                            |
| value_paths      | List[str]       | Paths to extract (supports * wildcard)              |
| exclude_patterns | List[str]       | Glob patterns to skip (e.g., ["EMG_*_metrics"])     |
| filters          | dict            | Optional filters from <code>create_filters()</code> |
| flatten_values   | bool            | Flatten nested dicts into columns                   |

```
df = extract_nested(
    profiles,
    base_path="sensor_metrics.emg",
```

```

    level_names=["date", "session", "side"],
    value_paths=["EMG_intensity.*", "EMG_rest_recovery.*"],
    exclude_patterns=["EMG_daily_metrics", "EMG_weekly_metrics"],
)
# Output: subject_id | date | session | side | EMG_intensity.mean_percent_mvc | ...

```

`extract_flat(profiles, base_path, filters=None, max_depth=10) -> pd.DataFrame` Flatten everything under a path into wide-format (one row per subject).

```

df = extract_flat(profiles, base_path="sensor_metrics.emg.EMG_weekly_metrics")
# Columns: subject_id, left.EMG_apdf.active.p10, left.EMG_apdf.active.p50, ...

```

---

## Path Resolution Functions

`resolve_path(data, path, default=None)` Get value from nested dict using dot-notation.

```

value = resolve_path(profile, "sensor_metrics.emg.EMG_weekly_metrics.left.EMG_apdf.active.p50")
# Returns: 12.5 (or None if path doesn't exist)

```

`path_exists(data, path) -> bool` Check if a path exists.

`list_keys_at_path(data, path) -> List[str]` List all keys at a given path.

```

keys = list_keys_at_path(profile, "sensor_metrics.emg")
# ['23-09-2025', '24-09-2025', 'EMG_weekly_metrics']

```

---

## Filtering

`create_filters(...) -> Dict[str, Any]` Create a filters dictionary for controlling extraction.

| Parameter                     | Type                         | Description                            |
|-------------------------------|------------------------------|--|
| <code>subject_ids</code>      | <code>List[str]</code>       | Include only these subjects            |
| <code>exclude_subjects</code> | <code>List[str]</code>       | Exclude these subjects                 |
| <code>groups</code>           | <code>List[str]</code>       | Filter by <code>meta_data.group</code> |
| <code>date_range</code>       | <code>Tuple[str, str]</code> | (start, end) in YYYY-MM-DD format      |
| <code>require_keys</code>     | <code>List[str]</code>       | Only include subjects with these paths |
| <code>custom_filter</code>    | <code>Callable</code>        | Custom (subject_id, profile) -> bool   |

```

from oh_parser import create_filters

filters = create_filters(
    subject_ids=["103", "104", "105"],
    date_range=("2025-09-01", "2025-09-30"),
    require_keys=["sensor_metrics.emg.EMG_weekly_metrics"],
)

df = extract_nested(profiles, ..., filters=filters)

```

## Usage Examples

### Example 1: Weekly EMG Summary (Wide Format)

```
from oh_parser import load_profiles, extract

profiles = load_profiles("/path/to/OH_profiles/")

df = extract(profiles, paths={
    # Left side
    "p10_L": "sensor_metrics.emg.EMG_weekly_metrics.left.EMG_apdf.active.p10",
    "p50_L": "sensor_metrics.emg.EMG_weekly_metrics.left.EMG_apdf.active.p50",
    "p90_L": "sensor_metrics.emg.EMG_weekly_metrics.left.EMG_apdf.active.p90",
    "rest_pct_L": "sensor_metrics.emg.EMG_weekly_metrics.left.EMG_rest_recovery.rest_percent",
    # Right side
    "p10_R": "sensor_metrics.emg.EMG_weekly_metrics.right.EMG_apdf.active.p10",
    "p50_R": "sensor_metrics.emg.EMG_weekly_metrics.right.EMG_apdf.active.p50",
    "p90_R": "sensor_metrics.emg.EMG_weekly_metrics.right.EMG_apdf.active.p90",
    "rest_pct_R": "sensor_metrics.emg.EMG_weekly_metrics.right.EMG_rest_recovery.rest_percent",
})
```

### Example 2: All Session-Level Data (Long Format)

```
from oh_parser import load_profiles, extract_nested

profiles = load_profiles("/path/to/OH_profiles/")

df = extract_nested(
    profiles,
    base_path="sensor_metrics.emg",
    level_names=["date", "session", "side"],
    value_paths=[
        "EMG_session.*",
        "EMG_intensity.*",
        "EMG_apdf.active.*",
        "EMG_rest_recovery.*",
    ],
    exclude_patterns=["EMG_daily_metrics", "EMG_weekly_metrics"],
)
# Columns: subject_id, date, session, side, EMG_session.duration_s, ...
```

### Example 3: Filter by Subjects and Date Range

```
from oh_parser import load_profiles, extract_nested, create_filters

profiles = load_profiles("/path/to/OH_profiles/")

filters = create_filters(
    subject_ids=["80", "81", "82"],
    date_range=("2025-09-23", "2025-09-26"),
)

df = extract_nested(
    profiles,
    base_path="sensor_metrics.emg",
    level_names=["date", "session", "side"],
```

```

    value_paths=["EMG_intensity.mean_percent_mvc"],
    exclude_patterns=["EMG_daily_metrics", "EMG_weekly_metrics"],
    filters=filters,
)

```

#### Example 4: Check Data Availability

```

from oh_parser import load_profiles, summarize_profiles, create_filters

profiles = load_profiles("/path/to/OH_profiles/")

# Check which subjects have EMG data
summary = summarize_profiles(profiles)
print(summary[summary['has_emg'] == True])

# Filter to only subjects with weekly EMG
filters = create_filters(
    require_keys=["sensor_metrics.emg.EMG_weekly_metrics.left"],
)
df = extract(profiles, paths={...}, filters=filters)

```

#### Example 5: Manual Path Navigation

```

from oh_parser import load_profiles, get_profile, resolve_path, list_keys_at_path, path_exists

profiles = load_profiles("/path/to/OH_profiles/")
profile = get_profile(profiles, "103")

# Check if path exists
if path_exists(profile, "sensor_metrics.emg.EMG_weekly_metrics"):
    print("Has weekly EMG!")

# Get a specific value
p50 = resolve_path(profile, "sensor_metrics.emg.EMG_weekly_metrics.left.EMG_apdf.active.p50")

# List available dates
dates = list_keys_at_path(profile, "sensor_metrics.emg")

```

---

## Path Syntax & Wildcards

### Dot Notation

Navigate nested structures with dots:

```
sensor_metrics.emg.EMG_weekly_metrics.left.EMG_apdf.active.p50
```

### Wildcards in value\_paths

Use `.*` to extract all keys under a path:

```

value_paths=["EMG_intensity.*"]
# Extracts: EMG_intensity.mean_percent_mvc, EMG_intensity.max_percent_mvc, ...

value_paths=["EMG_apdf.active.*"]
# Extracts: EMG_apdf.active.p10, EMG_apdf.active.p50, EMG_apdf.active.p90

```

## Level Names (Implicit Wildcards)

The `level_names` parameter creates wildcards for dynamic keys:

```
base_path="sensor_metrics.emg"
level_names=["date", "session", "side"]
# Iterates: sensor_metrics.emg.{date}.{session}.{side}
```

## Exclusion Patterns

Glob-style patterns to skip keys:

```
exclude_patterns=["EMG*_metrics"] # Skips EMG_daily_metrics, EMG_weekly_metrics
exclude_patterns=["*_aggregate"]  # Skips anything ending in _aggregate
```

---

## CLI Interface

Quick exploration without writing code:

```
# Basic info (default path)
python -m oh_parser --dir /path/to/OH_profiles

# List all subject IDs
python -m oh_parser --dir /path/to/OH_profiles --list

# Inspect a subject's profile structure
python -m oh_parser --dir /path/to/OH_profiles --inspect 103 --depth 5

# List all available paths
python -m oh_parser --dir /path/to/OH_profiles --paths 103

# Show data availability summary
python -m oh_parser --dir /path/to/OH_profiles --summary

# Quiet mode (suppress loading messages)
python -m oh_parser --dir /path/to/OH_profiles --quiet
```

---

## EMG Data Reference

### Metric Groups

#### EMG\_session - Session Metadata

| Key                            | Type  | Unit    | Description                      |
|--------------------------------|-------|---------|----------------------------------|
| <code>duration_s</code>        | float | seconds | Total recording duration         |
| <code>mvc_peak</code>          | float | mV      | MVC value used for normalization |
| <code>active_duration_s</code> | float | seconds | Time above rest threshold        |

#### EMG\_intensity - Intensity Metrics

| Key                  | Type  | Unit   | Description                       |
|----------------------|-------|--------|-----------------------------------|
| mean_percent_mvc     | float | %MVC   | Mean amplitude                    |
| max_percent_mvc      | float | %MVC   | Peak amplitude                    |
| min_percent_mvc      | float | %MVC   | Minimum amplitude                 |
| iemg_percent_seconds | float | %MVC*s | Integrated EMG (area under curve) |

#### EMG\_apdf - Amplitude Probability Distribution Function

| Key        | Type  | Unit | Description                               |
|------------|-------|------|---|
| full.p10   | float | %MVC | 10th percentile (all samples)             |
| full.p50   | float | %MVC | 50th percentile / median                  |
| full.p90   | float | %MVC | 90th percentile                           |
| active.p10 | float | %MVC | 10th percentile (active only, >=0.5% MVC) |
| active.p50 | float | %MVC | 50th percentile (active only)             |
| active.p90 | float | %MVC | 90th percentile (active only)             |

#### EMG\_rest\_recovery - Rest/Recovery Metrics

| Key                      | Type  | Unit     | Description                      |
|--------------------------|-------|----------|----------------------------------|
| rest_percent             | float | %        | Time below 0.5% MVC threshold    |
| gap_frequency_per_minute | float | gaps/min | Micro-break frequency            |
| max_sustained_activity_s | float | seconds  | Longest continuous active period |
| gap_count                | int   | count    | Total number of rest gaps        |


#### EMG\_relative\_bins - Relative Intensity Distribution

| Key              | Type  | Unit | Description           |
|------------------|-------|------|-----------------------|
| below_usual_pct  | float | %    | Time below weekly P10 |
| typical_low_pct  | float | %    | Time between P10-P50  |
| typical_high_pct | float | %    | Time between P50-P90  |
| high_for_you_pct | float | %    | Time above weekly P90 |


#### Aggregation Levels

| Level   | Path                                | Description                       |
|---------|-------------------------------------|-----------------------------------|
| Session | emg.{date}.{session}.{BinId}        | Recording metrics                 |
| Daily   | emg.{date}.EMG_daily_metrics{BinId} | Durations-weighted daily average  |
| Weekly  | emg.EMG_weekly_metrics{BinId}       | Durations-weighted weekly average |

#### Pipeline Architecture



## Data Flow



## Notes

- **Date format:** EMG paths use DD-MM-YYYY format
- **Time format:** Session times use HH-MM-SS format
- **Null values:** Some metrics may be **None** (e.g., relative bins at weekly level)
- **Sides:** EMG data has separate entries for **left** and **right**
- **macOS:** Hidden **.\_\*** files are automatically skipped when loading

---

*Generated for oh\_parser module - PrevOccupAI+ Project*